

Types of Hearing Aids

Hearing aids have come a long way, offering a variety of styles to match every lifestyle and hearing need. Whether you're an active adventurer, a busy professional, or someone who enjoys quiet moments at home, there's a hearing aid that's perfect for you.



Completely-in-Canal (CIC)

The most discreet option, CIC hearing aids sit deep within the ear canal. They are perfect for people who want their hearing aids to be completely unnoticeable.



In-the-Canal (ITC)

These tiny devices fit partially in the ear canal, offering a nearly invisible solution. They are custom-made to best fit you and are suited for individuals who prioritise aesthetics and convenience.



In-the-Ear (ITE)

These custom-made devices fit in the ear and are best for individuals who prioritise aesthetics and convenience.



Receiver-in-Canal (Slim - RIC)

Hearing aids that use a thin wire that extends from the hearing aid, over the outer ear, and into the ear canal for a natural sound experience.



Receiver-in-Canal (RIC)

Hearing aids that use a thin wire that extends from the hearing aid, over the outer ear, and into the ear canal for a natural sound experience.



Behind-the-Ear (BTE)

Hearing aids that sit behind-the-ear with the impressive aesthetics, comfortable fit and intuitive controls. These classic models are known for their durability and power.